

Lubtear[®] Eye Drops

Dextran 70 USP 0.1% & Hypromellose USP 0.3% Sterile Eye Drops

Dextran 70 USP & Hypromellose USP combination Eye drops acts like natural tears, which provides smoothing relief to dry, itchy, burning and irritated eyes. It helps eyes to stay moist healthy and protects against further irritation. This unique combination of ingredients works together to retain moisture on the eye and slow evaporation of the tear film.

COMPOSITION

Lubtear[®] Eye Drops: Each ml contains Dextran 70 USP 1 mg and Hypromellose USP 3 mg.

USES & INDICATION

As a lubricant and artificial tear in dry eye and other ocular irritation syndromes associated with deficient tear or mucous secretion. This combination also prevents cornea to damage in patients with keratoconjunctivitis and for ocular lubrication. It is also used for the temporary relief of burning and irritation due to dryness of the eye and for use as the protectant against further irritation.

DOSAGE & ADMINISTRATION

Adults and children: One or two drops three times daily or as directed by the physician.

CONTRAINDICATION

This product contains Benzalkonium Chloride BP and should not be used when soft contact lenses are being worn.

PRECAUTION

If irritation persists or if the condition does not improve, patient should seek further advice from doctor.

USE IN PREGNANCY & LACTATION

There is insufficient evidence as to the safety in pregnancy and lactation. Therefore, this product should only be used in pregnancy and lactation if it is considered essential by the physician.

SIDE EFFECT

There are no known side effects with the use of it, however, if the patient experiences any reaction in eye or other part of the body after using this medication then consult with doctor.

STORAGE

Store below 30°C and do not freeze. Do not use after 28 days of first opening. Keep the container closed tightly after each opening. Keep out of reach of children.

HOW SUPPLIED

Lubtear[®] Eye Drops: Each plastic dropper bottle contains 10 ml of Dextran 70 USP & Hypromellose USP combination sterile solution.

SQUARE